
VenaTech® LP

REF. 04435125

F

B | BRAUN

CE
0459

ENGLISH en

FRANÇAIS fr

DEUTSCH de

ESPAÑOL es

ITALIANO it

SVENSKA sv

PORTUGUÊS pt

NORSK no

DANSK da

SUOMI fi

POLSKI pl

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
VenaTech[®] LP

REF. 04435125

Instructions for use for implantation via the :
femoral approach


**Use this booklet only if you have chosen to implant the filter
from the femoral approach.**

Do not use this booklet for any other approach.

B | BRAUN

System contents

- (1) Filter in its cartridge
- (2) Cartridge cradle
- (3) Introducer sheath with hemostasis valve and side connection
- (4) Dilator
- (5) Long Pusher (white)
- (6) Short Pusher (blue)
- (7) PTFE coated "J" guide, length = 150cm, diameter = 0.89mm (0.035 in)

Important warnings

- **Instruction manual**

This instruction manual is to be used **when implantation from the femoral approach** has been selected. Set aside the manual provided for the other approaches prior to beginning the implantation procedure.

- **Implantation period**

This filter is designed to provide effective protection against **permanent or long term risk of pulmonary embolism**.

- **Implantation approach**

This filter is designed for insertion via :

- the right femoral vein
- the left femoral vein

when placement from below has been selected.

Warning : Avoid the use of a venous access site which was previously used for an implanted central venous catheter. Implantation of a vena cava filter using an existing access site can result in incomplete filter deployment. This could in turn result in filter migration and/or inadequate protection against pulmonary embolism.

Note : Extreme care must be exercised during approaches through tortuous anatomies which can cause sheath kinking and make filter insertion difficult or impossible. Under extreme conditions, attempting to force the filter through a kinked sheath could result in a sheath perforation. In any cases, stop pushing on the filter when encountering significant resistance.

- **Diameter of the vena cava**

The maximum diameter of the vena cava, as evaluated during cavography, must be less than or equal to **35 mm** from an AP projection when corrected for magnification.

In cases in which measurement is accomplished using Computed Tomography (CT), the section with the largest diameter should be used.

Important : Do not attempt to increase the diameter of the filter.

- **Filter placement in pregnant women**

In pregnant patients where the fetus may be endangered by fluoroscopy, risks and benefits should be carefully assessed.

• Compatibility with Magnetic Resonance Imaging (MRI)

Conditional

The VenaTech® LP, Vena Cava Filter was determined to be MR-Conditional according to the terminology specified in the American Society for Testing and Materials (ASTM) International, Designation: F2503-08. Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment.

Non-clinical testing demonstrated that the VenaTech® LP, Vena Cava Filter is MR Conditional. A patient with this device can be scanned safely, immediately after placement under the following conditions:

- Static magnetic field of 3-Tesla and 1.5-Tesla
- Maximum spatial gradient magnetic field of 710 Gauss/cm or less
- The scanner can be operated in the Normal Operating Mode with a maximum whole body averaged specific absorption rate (SAR) of 2.0 W/kg or in the First level Controlled Mode with a maximum SAR of 4.0 W/kg for 15 minutes of scanning.

MRI-Related Heating

In non-clinical testing, the VenaTech® LP Vena Cava Filter produced the following temperature rise during MRI performed for 15-min (i.e., per pulse sequence) in 3-Tesla (Excite, Software G3.0-052B, General Electric Healthcare, Milwaukee, WI) MR systems, as:

MRI Condition	MR System Reported, Maximum Whole Body Average SAR (W/kg)	Calorimetry Value (W/kg)	Highest Temperature Change	Time for MRI (per pulse sequence)
3-T / 128-MHz	2.9	2.7	+ 2.8 °C	15-min

Artifact Information

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the VenaTech® LP, Vena Cava Filter. Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary.

The maximum artifact size extends approximately 11-mm relative to the size and shape of the VenaTech® LP, Vena Cava Filter when scanned in non-clinical testing using the sequence gradient echo pulse sequence in a 3T MR system (Intera Achieva upgrade) – Phillips Medical Systems Best, The Netherlands (software version: 2.5.3.3 2007-09-28) with a transmit/receive birdcage coil, 28 cm long, 27 cm in diameter.

- **Secondary Catheterization**

Following filter placement, no central venous catheterization should be attempted without fluoroscopic guidance.

Potential Adverse Effects

- Thrombosis at puncture site
- Perforation of the inferior vena cava (IVC) and adjacent organs
- Migration (localized or intra-cardiac)
- Filter breaking
- Filter tilting
- Recurrence of pulmonary embolism under efficient filtration
- Thrombosis of the IVC
- Venous insufficiency
- Death

Indications

The VenaTech® LP filter is designed for patients in whom the risk of thromboembolism is believed to be permanent or of very long duration.

Indications for filter implantation are the following :

- Patients for whom anticoagulants are contraindicated
- Patients for whom anticoagulation resulted in complications
- Patients for whom anticoagulants failed or were insufficient
- Emergency treatment following massive pulmonary embolism where anticipated benefits of conventional therapy are reduced

Contraindications

- Vena cava filters should not be implanted in patients with **risk of septic embolism**.
- Patients with a vena cava which has a **diameter greater than 35 mm** (due to the risk of device migration).
- Known allergy to the materials contained in the kit (main raw materials: $\text{CoCr}_{20}\text{Ni}_{16}\text{Mo}_7$, Polycarbonate, Polysulfone, Silicone, Tetrahydrofurane, Methylene Chloride, Polyvinyl chloride, Polyethylene, Fluorinated Ethylene Propylene, Stainless steel, Polytetrafluoroethylene).

Implantation procedure

Pre-implant cavography is required :


- To confirm the patency and visualize the anatomy of the vena cava,
- To mark the level of the renal veins,
- To locate the highest level of any thrombus which may be present,
- To determine the desired level for filter deployment and to mark the position with respect to the vertebral bodies,
- To confirm that the diameter of the vena cava (AP projection) at the site where the filter is to be deployed is less than or equal to the maximum authorized diameter.

A) Preparation of the implantation accessories

1. The filter is pre-loaded into its cartridge. Load the cartridge into the cradle as shown on the photo.

Warning 1 : A “click” must be heard to guarantee complete insertion.

Warning 2 : The side window of the cradle must read “FEMORAL”.

2. Immerse the assembly in a container of **heparinized** saline.

Note : This precaution reduces the possibility of thrombus formation around the filter during its movement through the introducer sheath.

B) Puncture of the vein

1. Prepare the patient for access from the selected approach.
2. Puncture the vein using the Seldinger technique (or perform a surgical cut-down of the vein). Confirm venous reflux into the syringe.
3. Advance a guidewire (0,89 mm – 0.035 in.) and then remove the Seldinger needle.

C) Control cavogram, marking the filter implant site


1. Create a small incision with a scalpel to facilitate the passage of the introducer. Insert a 7F catheter introducer over the guidewire. Perform a cavogram using the 7F catheter introducer or a catheter of your choice.

Caution : If ilio-caval clot is present, extreme care must be exercised to avoid dislodging the clot.

2. Locate the renal veins and select the location for filter placement. Mark the planned position of the top of the filter cone by placing a metallic marker (e.g. a hemostat) horizontally across the abdomen.

Note : Choice of the site for filter deployment :

When possible, it is preferable to place the filter just below the renal veins. A security margin of 0.5 cm should be left from the level of the lowest renal vein to the top of the filter.

D) Passage of the "J" guidewire

Advance the "J" guidewire provided into the kit (150 cm - diameter 0,89 mm / 0.035 in.). Follow its progress under fluoroscopy.

Warning : To prevent the introducer from inadvertently passing outside of the IVC (e.g. into a tributary vessel) verify that the "J" portion of the guidewire is located above the renal veins.

Note : In the event of difficulty advancing the guidewire because of the "J", remove the guidewire and use the other end, which is straight and supple.

E) Placement of the introducer


1. Insert the dilator into the introducer sheath. Connect the Luer-Lock connecting the sheath and dilator.
2. Place the introducer assembly over the "J" guidewire. Advance it over the guidewire.
3. Continue to advance the introducer until the marker is a little bit further than one cm above the deployment site marker (hemostat) placed on the abdomen.

Note : The level of the marker on the introducer corresponds to the top of the filter cone when it is collapsed within the sheath prior to deployment.

This maneuver pre-positions the filter slightly above the level chosen for deployment. Withdrawing the sheath slightly adjusts the filter position before deployment.


Note : Final adjustment of the filter position will be made just prior to deployment.

Warning : Hold the guidewire firmly as the introducer is advanced to prevent the guidewire from moving or developing a kink.

F) Placement of the filter into the sheath

1. Disconnect the luer lock connecting the sheath and dilator.
2. Remove the "J" guidewire and the dilator taking care not to move the sheath.
3. Flush through the side connection of the sheath with 10 ml of heparinized saline.
4. Screw the cartridge cradle with its filter clockwise onto the sheath. Ascertain that both are securely locked together.
5. Using the short pusher (blue), insert the filter into the sheath. The pusher must be completely inserted.
6. Unscrew the cartridge cradle and remove it with the pusher as a unit.

Note : The cartridge cradle must be removed to ensure a correct filter implantation.

7. Inject several milliliters of heparinized saline to prevent thrombus formation in the sheath, using the side connection.

G) Advance the filter through the sheath

1. Insert the distal end of the long pusher (white) into the sheath. Hold the sheath firmly, not allowing it to move.
2. Using the long pusher (white), advance the filter through the sheath under fluoroscopy, until the proximal mark on the pusher is aligned with the mouth of the introducer sheath.
3. The filter is now located at the distal end of the sheath. The system is ready for final positioning prior to deployment.



H) Deployment of the filter

1. **Under fluoroscopic guidance, align the top of the filter cone with the radiopaque marker** placed on the patient at the beginning of the procedure. (This is the final adjustment of filter position prior to deployment).
2. Once the filter has reached the proper level for deployment, **be careful not to move the introduction system.**
3. Confirm under fluoroscopy that the sheath has not moved and that the filter is properly positioned at the desired level for deployment. The top of the filter should be at the level of the metallic marker placed on the patient.
4. Holding the pusher firmly with one hand, take the sheath in the other hand and, under fluoroscopic guidance, withdraw the sheath. This movement deploys the filter.

Reminder : To assure good filter position, it is extremely important not to move the pusher :

- during the step prior to filter deployment,
- during filter deployment.

5. The filter is now implanted.

Note : To perform a cavogram at this time,

1. Remove the pusher while keeping the sheath in place.
 2. Perform the cavogram by manually injecting contrast material through the side connection.
6. Remove the assembly sheath/pusher (or sheath alone if the pusher has already been removed to perform a cavogram).

I) Final check

Take a final film which documents the location and deployment of the filter.

Important

Vena cava filter for permanent interruption of the Inferior Vena Cava, STERILE and non-pyrogenic, sterilized using Ethylene Oxide. Check that the sterile packaging is undamaged before opening. For single use only. The device and its accessories are not re usable nor designed to be re-used. Any re-use would definitely compromise the performance and safety of the device. Destroy system components after filter implantation. Components may not be cleaned, resterilized or reused. Store at room temperature. Avoid freezing and excessive heat. Do not remove from unit carton until ready for use.

Please review the sections of the instructions for use which describe :

- the implantation procedure
- the indications for filter placement
- the contraindications to filter placement.

Opening the sterile packaging :

Confirm the date of sterile expiration on the labelling. Never implant an IVC filter which has exceeded its sterile expiration date. To maintain sterility, all normal operating room sterile procedures should be followed.

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